

## Product Datasheet

Model: T14M

Description:



14x17" MOBILE TRANS-FORMER™ DRP HOLDER

Revision: MARCH 2014

### Product Attributes:

- Holds DR Panel during exams.
- Mobile.
- Protects grid from incidental damage during use and storage.
- Provides additional protection for DR Panel.
- Cabinet adjusts vertically, and pivots 90 degrees in either direction.
- Features removable grid.

### Technical Specifications

Model	T14M
Weight (approx.)	86 lbs (39kg)
Size	22.5" x 22.5" x 70" (57.1cm x 57.1cm x 177.8cm)
Vertical Travel (height adjustment)	Top - 75.75" (192.4cm)
Boxed Dimension	31"x31"x80" (78.74cm x 78.74cm x 203.2 cm)
Shipping Weight	130 lbs. (59 kg)

### Weight Tolerance

*Note: This product was not designed for weight-bearing studies and should never be used without the DR Panel installed.*

Weight Tolerance: 400 lbs. (182 kg)

### Flammability Rating

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

### Biocompatibility

This model is a biocompatible device per ISO 10993.

### Statement of Conformance

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

### Registrations and Compliance

All models listed above are made in the USA and registered to:

**Reina Imaging**

**6107 Lou Ave.**

**Crystal Lake, IL 60014 USA**

### Device Family Classification

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE

TGA	Class 1
-----	---------