

Mobile DR Panel Full Enclosure

Product Advantages

IP 5X Rating

Shock absorbent foam lining

Quick release locks

Superior particulate and dust protection

Smooth non-porous surface

Curved corners

Available in 14x17in and 17x17in



The **14x17 Mobile DR Panel Full Enclosure** provides additional fluid and drop protection to sensitive flat panel detectors. With **quick release locks**, the case is both efficient and easy to use allowing the radiographer to focus on **positioning and patient care** instead of panel protection. Ideal for the ER, OR, disaster relief and any environment that is subject to potential fluid damage or accidental drops.

ReinaImaging | XRAY CASSETTE REPAIR

Reimagine Radiography



www.reinaimaging.com

800-752-4918

customer_service@reinaimaging.com

Product Details



Tube Side



Rear Side



Cutouts for LEDs and PB



Quick Release Locks

Technical Data for FE14

Product Dimensions

19.34 x 16.71 x 0.944 in. (49.12 x 42.44 x 2.39 cm)

Weight

3.20 lbs. (1.45 kg)