

PRODUCT DATASHEET



Model: FE14

Description: **14x17" Mobile DR Panel Full-Enclosure**

Revision: March 2017

PRODUCT ATTRIBUTES:

- Fully encases DR Panel in protective plastic shell.
- DR Panel is removable for charging.
- Designed for and fits most wireless 14x17" cassette sized DR Panels.
- Portable.



MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Shell	HDPE and Polycarbonate	98
Overlay	Polycarbonate Film	2

DIMENSIONS AND PERFORMANCE

Product Dimension: 19.34 x 16.71 x 0.944 in. (49.12 x 42.44 x 2.39 cm)

Weight: 3.20 lbs. (1.45 kg)

Weight Tolerance: 400 lbs. (181.4 kg)

Not designed for weight bearing exams. Not to be used without an image receptor.

IP RATING

IP5X for dust per IEC 60529, Ed. 2.2, Sec. 13.4

FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging
6107 Lou Ave.
Crystal Lake, IL 60014 USA

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.
20148Milano Via don Gnocchi 24, Italy
+39 02 48706799

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

6107 West Lou Avenue | Crystal Lake, IL 60014 | (P) 800.752.4918 | (F) 815.356.8270 | reinaimaging.com