

PRODUCT DATASHEET

Model: CC5NG17, CC5NG14, CC5NG10, CC5NG24

Description: CarbonLite Mobile DR Panel Holder

Revision: June 2020



Product Attributes

- CarbonLite™ features .5mm carbon fiber
- Single ergonomic handle on long axis
- Rigid glass-filled polypropylene frame
- Curved corners and radiused edges
- Smooth, non-porous surfaces
- Simple Twin-Lock™ dual locking design
- Shock absorbent corner bumpers and inserts
- Fits all 17x17in, 14x17in, 10x12in and 24x30cm ANSI standard sizes

Material Composition

Material	%
Frame	
GF Polypropylene	85%
Aluminum	6%
Stainless Steel	5%
Flex PVC	2%
Acetal	1%
PC ABS	1%
Graphics Overlay	
Polycarbonate	98%
Colorant	2%
Tubeside Protection	
Carbon Fiber	98%
Polycarbonate	2%

PRODUCT DATASHEET

Product Code	Product Name	Dimensions (cm)	Weight (kg)
17 x 17in			
CC5NG17	17x17 CarbonLite Mobile DR Panel Holder	56.42 x 49.22 x 2.69	1.76
14 x 17in			
CC5NG14	14x17 CarbonLite Mobile DR Panel Holder	48.91 x 49.22 x 2.69	1.58
10 x 12in			
CC5NG10	10x12 CarbonLite Mobile DR Panel Holder	38.47 x 36.54 x 2.69	1
24 x 30cm			
CC5NG24	24x30 CarbonLite Mobile DR Panel Holder	38.47 x 36.54 x 2.69	1

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein in this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:
Reina Imaging
6107 Lou Ave.
Crystal Lake, IL 60014 USA