

PRODUCT DATASHEET

Model: CC2NG17, CC2NG14, CC2NG10, CC2NG24

Description: CarbonMax Mobile DR Panel Holder

Revision: June 2020



Product Attributes

- CarbonMax™ features 2.5mm carbon fiber
- Single ergonomic handle on long axis
- Rigid glass-filled polypropylene frame
- Curved corners and radiused edges
- Smooth, non-porous surfaces
- Simple Twin-Lock™ dual locking design
- Shock absorbent corner bumpers and inserts
- Fits all 17x17in, 14x17in, 10x12in and 24x30cm ANSI standard sizes

Material Composition

| Material | % |
|-----------------------------|------|
| Frame | |
| GF Polypropylene | 85% |
| Aluminum | 6% |
| Stainless Steel | 5% |
| Flex PVC | 2% |
| Acetal | 1% |
| PC ABS | 1% |
| Graphics Overlay | |
| Polycarbonate | 98% |
| Colorant | 2% |
| Tubeside Protection | |
| Carbon Fiber | 100% |

PRODUCT DATASHEET

| Product Code | Product Name | Dimensions (cm) | Weight (kg) |
|------------------|--|----------------------|-------------|
| 17 x 17in | | | |
| CC2NG17 | 17x17 CarbonMax Mobile DR Panel Holder | 56.42 x 49.22 x 2.69 | 1.76 |
| 14 x 17in | | | |
| CC2NG14 | 14x17 CarbonMax Mobile DR Panel Holder | 48.91 x 49.22 x 2.69 | 1.54 |
| 10 x 12in | | | |
| CC2NG10 | 10x12 CarbonMax Mobile DR Panel Holder | 38.47 x 36.54 x 2.69 | 1 |
| 24 x 30cm | | | |
| CC2NG24 | 24x30 CarbonMax Mobile DR Panel Holder | 38.47 x 36.54 x 2.69 | 1 |

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA