

PRODUCT DATASHEET

Model: FS

Description: **WEIGHT-BEARING FOOT STOOL**

Revision: May 2017

PRODUCT ATTRIBUTES:

- Lightweight.
- Handrail.
- Slip resistant base.
- Accommodates 14x17" cassette and smaller.
- Low to ground.

DIMENSIONS AND PERFORMANCE

Base Dimensions: 22 x 22 x 9.25 in. (55.9 x 55.9 x 23.5 cm)

Handle Height to Floor: 45 in. (114.3 cm)

Weight: 33.6 lbs. (15.24 kg)

Weight Tolerance: 500 lbs. (226.8 kg)

Cassette/Panel Slot Dimensions: 22.25 x 0.75 x 2.125 in. (56.5 x 1.9 x 5.4 cm)

Boxed Dimensions and Weight

Box 1: 24x24x12 in. (61x61x30.5 cm), 35 lbs. (15.87 kg)

Box 2: 48x20x3 in. (122x50.8x7.6 cm), 11 lbs. (4.99 kg)



FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein in this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.


Crystal Lake, IL 60014 USA

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148 Milano Via don Gnocchi 24, Italy

+39 02 48706799



DEVICE FAMILY CLASSIFICATION

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1