

PRODUCT DATASHEET

Model: P2C10

Description: Protect-A-Grid 2 With Channels 10x12"

Revision: July 2020



Product Attributes

- Lightweight
- Thin profile
- Curved corners
- Economical
- Works with all 10x12in. film and CR cassettes

Material Composition

Material	%
Shell, Aluminum	98%
Overlay, Polycarbonate Film	2%

DIMENSIONS

Product Dimension: 11.38x13.06x0.88in. (28.89x33.18x2.22cm)

Weight: 1.65lb. (0.75kg)

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA