

# PRODUCT DATASHEET

Model: PB3HNG14

Description: **HD Mobile DR Panel Holders 3-Handle**

Revision: September 2019



## MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Shell	HDPE and Polycarbonate	98
Overlay	Polycarbonate Film	2

## DIMENSIONS AND PERFORMANCE

Product Dimension: 51.4 x 45 x 2.54 cm. (20.25 x 17.75 x 1.0 in)

Weight: 1.34 kg (2.751lb.)

Weight Tolerance: 400 lbs. (181.4 kg)

*Not designed for weight bearing exams. Not to be used without an image receptor.*

## FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

## BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

## STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

## REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

**Reina Imaging**

6107 Lou Ave.

Crystal Lake, IL 60014 USA

## AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148 Milano Via don Gnocchi 24, Italy

+39 02 48706799

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1