

Product Datasheet

Model: PS1417
 Description: PROTECT-A-GRID SNAP-ON ENCASUREMENT
 Revision: MAY 2013



Product Attributes:

- Encases grid to protect against damage.
- Holds DR Panel during exams using locking latches.
- Provides additional protection for DR Panel.

Dimensions

Product Dimension: 19 x 15.5 x 1" (48.26 x 39.37 x 2.54 cm)
 Weight: 0.88 lbs. (0.4 kg)*

Material Composition

Part	Material	Approximate %
Shell	HDP	98
Overlay	Polycarbonate Film	2

Weight Tolerance

Note: This product was not designed for weight-bearing studies and should never be used without the DR Panel installed.

Weight Tolerance: 400 lbs. (181.4 kg)

Flammability Rating

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

Biocompatibility

This model is a biocompatible device per ISO 10993.

Statement of Conformance

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

Registrations and Compliance

All models listed above are made in the USA and registered to:

Reina Imaging
 6107 Lou Ave.
 Crystal Lake, IL 60014 USA

Authorized EU Representative

PINCHET RX. S.r.l.
 20148 MILANO, Via don Gnocchi 24, ITALY

Device Family Classification

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

*Weight without installed grid.