

# PRODUCT DATASHEET

Model: **WDDR14-H**

Description: **WEIGHT DISTRIBUTION CASSETTE CAP WITH HANDLE FOR 14x17 (35x43cm) DR PANELS**

Revision: May 2019

## PRODUCT ATTRIBUTES:

- Protects DR Panel during weight-bearing exams.
- Fits over the DR Panel.
- Features slip resistant grips for stability.
- For non-gridded weight-bearing studies.
- Portable. Features a stainless steel handle.



## DIMENSIONS AND PERFORMANCE

Product Dimension: 21.5 x 16.875 x 2 in. ( 54.61 x 42.86 x 5.08 cm)

Weight: 7.2 lbs. (3.26 kg)

Weight Tolerance: 750 lbs. (340 kg) with the Detector Installed

## FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel. UL Rating: UL94V-2.

## BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

## STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

## REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

**Reina Imaging**

**6107 Lou Ave.**

**Crystal Lake, IL 60014 USA**

## AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

DEVICE FAMILY CLASSIFICATION	
Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

MATERIAL COMPOSITION		
PART	MATERIAL	APPROXIMATE %
Shell	HDP	95
Overlay	Polycarbonate Film	4
Handle	Stainless Steel	1