

PRODUCT DATASHEET

Model: **WDDR17-C**


Description: **WEIGHT DISTRIBUTION CASSETTE CAP WITH CLEAR TOP FOR 17x17 (43x43cm) DR PANELS**

Revision: April 2018

ReinaImaging

PRODUCT ATTRIBUTES:

- Protects DR Panel during weight-bearing exams.
- Fits over the DR Panel with clear top.
- Features slip resistant grips for stability.
- For non-gridded weight-bearing studies.
- Portable.



MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Shell	Polycarbonate and Polypropylene	98
Overlay	Polycarbonate Film	2

DIMENSIONS AND PERFORMANCE

Product Dimension: 20.125 x 20.125 x 2.25 in. (51.1 x 51.1 x 5.7 cm)

Weight: 8.4 lbs. (3.81 kg)

Weight Tolerance: 750 lbs. (340 kg)

FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein in this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

DEVICE FAMILY CLASSIFICATION

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

6107 West Lou Avenue | Crystal Lake, IL 60014 | (P) 800.752.4918 | (F) 815.356.8270 | reinaimaging.com