



Product Datasheet  
Models: T34-C, T34M-C  
Revision: Original (March 2011)

### Trans-Former Cabinet for use with the Carestream LLI Grid Assembly

Product Attributes:

- Mobile or Wall-Mounted adjustable cabinet for Kodak LLI Grid
- Adjusts vertically 35.5" (90 cm) from lowest to highest position
- Non-electrical
- Mobile model features casters
- Positive lock height adjustment

### Dimensions

Product Dimension	21.25" x 16.5" x 53.5" (53.975 cm x 41.91 cm x 135.89 cm) (Length x Width x Height)
Weight (without LLI Grid)	84 lbs (38.1 kg)
Boxed Dimensions	31" x 31" x 60" (78.74 cm x 78.74 cm x 152.4 cm) (Length x Width x Height)
Shipping Weight	120 lbs. (54.43 kg)



### Flammability Rating

The product is non-flammable under the normal operating conditions of the Carestream LLI Grid.  
UL Rating: UL94V-2

### Biocompatibility

The T34-C and T34M-C are biocompatible devices.

### Registrations and Compliance

Product registered to:

**Reina Imaging**  
**6107 Lou Ave.**  
**Crystal Lake, IL 60014**

- T34M-C and T34-C are Class 1 Medical Devices in the USA
- 510(k) exempt
- Made in the USA
- CE