

PRODUCT DATASHEET

Model: **T34M-DV, T34-DV**

Description: **TRANS-FORMER LLFS GRID CABINET FOR CARESTREAM DIRECTVIEW CR CASSETTES**

Revision: May 2018

PRODUCT ATTRIBUTES:

- Mobile or wall mounted
- Houses a removable grid cabinet for Directview CR Cassettes.
- Encases grid in protective plastic to protect against damage.
- Cabinet adjusts vertically and pivots 90 degrees.
- Non-electrical.

MOBILE DIMENSIONS AND PERFORMANCE

Product Dimension: 21.25" x 16.5" x 53.5" (53.975 cm x 41.91 cm x 135.89 cm)

Vertical Travel (Height Adjustment): 35.5" (90 cm)

Weight: 86 lbs. (39 kg)

Weight Tolerance: 400 lbs. (182 kg) with Cassettes Installed (not designed for weight bearing exams)

Boxed Dimension: 31"x31"x60" (78.74cm x 78.74cm x 152.4 cm)

Shipping Weight: 120 lbs. (54.43 kg)

WALL-MOUNTED DIMENSIONS AND PERFORMANCE

Product Dimension: 22.5 x 53.6 in. x 6" (57.15 x 136.14 cm x 15.24 cm)

Vertical Travel (Height Adjustment): Top – 74.75" (189.86 cm.)

Weight: 50 lbs. (22.67 kg)

Boxed Dimension: 60 x 26 x 12 in. (152.4 x 66.04 x 30.48 cm)

Shipping Weight: 70 lbs. (31.75 kg)

FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.
UL Rating: UL94V-2.

BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein in this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA



Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1