

Product Datasheet

Model: T50M-3C
Description: TRANS-FORMER™ MOBILE GRID CABINET FOR 3 DR PANEL LONG-LEG/FULL-SPINE IMAGING SYSTEMS
Revision: OCTOBER 2015

Product Attributes:

- Mobile.
- Houses a removable grid cabinet for 3 DR Panel Imaging Systems.
- Encases grid in protective polycarbonate plastic to protect against damage.
- Cabinet adjusts vertically.
- Non-electrical.

Technical Specifications

Model	T50M-3C
Weight (approx.)	86 lbs (39 kg)
Size	22.5 x 53.6 in. (57.15 x 136.14 cm)
Vertical Travel (Height Adjustment)	Top – 75.75" (192.4 cm.)
Boxed Dimension	31 x 31 x 60 in. (78.74 x 78.74 x 152.4 cm)
Shipping Weight	130 lbs. (59 kg)

Weight Tolerance

Note: This product was not designed for weight-bearing studies and should never be used without the CR Cassettes installed.

Grid Cabinet
Weight Tolerance: 600 lbs. (272 kg) with cassettes installed

Flammability Rating

The product is non-flammable under normal operating conditions.
UL Rating: UL94V-2.

Biocompatibility

This model is a biocompatible device per ISO 10993.

Statement of Conformance

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

Registrations and Compliance

All models listed above are made in the USA and registered to:

Reina Imaging
6107 Lou Ave.
Crystal Lake, IL 60014 USA

Authorized EU Representative

PINCHET RX. S.r.l.
20148 MILANO, Via don Gnocchi 24, ITALY

Device Family Classification

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1