

PRODUCT DATASHEET

Model: **PDR710**

Description: **Protect-A-Grid DR for Canon CXDI -710C**

Revision: February 2018

PRODUCT ATTRIBUTES:

- Encases grid to protect against damage.
- Protects DR Panel during routine handling.
- Provides ergonomic handle and double positive lock to secure DR Panel during transport.
- Made of shock absorbing materials.

ReinaImaging



MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Shell	HDP	98
Overlay	Polycarbonate Film	2

DIMENSIONS AND PERFORMANCE

Product Dimension: 18.380 x 18.630 x 0.790 in. (46.68 x 47.32 x 2.0 cm)

Weight: 2.0 lbs. (0.907kg)

Weight with 8:1/132 lpi grid installed: 4.35lbs. (1.97 kg)

Weight Tolerance: 400 lbs. (181.4 kg)

Not designed for weight bearing exams. Not to be used without an image receptor.

FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

DEVICE FAMILY CLASSIFICATION

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

6107 West Lou Avenue | Crystal Lake, IL 60014 | (P) 800.752.4918 | (F) 815.356.8270 | reinaimaging.com