

# Clear Top Mobile Table

## Product Advantages

Clear polycarbonate top with scratch-resistant overlamine film

4 durable swivel casters with foot locks

All steel structure and welded base provides maximum strength and reliability

500 lb. (227kg) patient weight limit




The **Clear Top Table** is a mobile radiographic table that is a perfect fit for a variety of x-ray equipment systems used in hospitals, imaging centers, ortho clinics and private practice settings. The table **pairs perfectly with our line of Positioning Partners** and an OTC or floor mounted tube stand as an economic imaging table. With its four 360 degree swivel locking casters, heavy duty steel frame and radiolucent scratch resistant polycarbonate top, this elegantly simple table provides a versatile, reliable and affordable solution.

# Product Details



\*Positioning Partner/Imaging receptor not included

# Product Details



## Technical Data

Table top length	80 in. (203.2cm)
Table top width	26 in. (66.1cm)
Table top height	27.9 in. (70.9cm)
Clear span, length	68 in. (172.7cm)
Clear span, width	19.25 in. (48.9cm)
Overall Weight	164 lbs. (74kg)
Weight Capacity	500 lb. (227kg)