



EASY READER™ RADIOPAQUE RULER

Product Advantages


- Effortless application
- Lead-free, environmentally friendly, safe for patient contact
- Stainless steel plastic-lined bracket slides easily across entire image area
- Clearly visible, high-contrast, easy to read radiopaque reference marks
- Made from resilient material that is both flexible and strong
- Available in various lengths with and without mounting bracket

The **Easy Reader** radiopaque ruler is ideal for spine and hips to ankle examinations. Its high-contrast NIST compliant reference marks are easy to see with your eyes as well as in a radiograph. It is simple to attach and remove with a stainless steel plastic-lined bracket that glides across the entire imaging area and won't scratch or mar equipment. Available in a variety of lengths, with or without a mounting bracket. It is 100% lead-free, RoHS compliant and made in the USA.





Hook allows for easy application, positioning and removal



Floor mounted column with long format detector



Attaches to Stitching Platform