

PRODUCT DATASHEET



Model: ER (see below for details)

Description: **Easy Reader Radiopaque Ruler**

Revision: August 2019

PRODUCT ATTRIBUTES:

- Hanging application
- Lead-free, environmentally friendly, safe for patient contact
- Clearly visible, high-contrast, easy to read radiopaque reference marks
- Made from resilient material that is both flexible and strong
- Available in various lengths with and without mounting bracket

MATERIAL COMPOSITION

| PART | MATERIAL | APPROXIMATE % |
|--------|----------|---------------|
| Shell | Lexan | 98 |
| Filler | Epoxy | 1 |
| Text | Tungsten | 1 |

DIMENSIONS AND PERFORMANCE

Product Dimension: See below

| Part | Description | Dimensions | Box Dim's & Weight (lb) |
|---------|---------------------------|--------------------------|-------------------------|
| ER66SPH | SPH | 1747 mm x 40 mm x 100 mm | 4x4x73" - 3.75 lb |
| ER34B | POSPIT2 | 1010 mm x 40 mm x 100 mm | 4x4x44" - 3.50 lb |
| ER34BT | 2 Panel Trans-Former | 1010 mm x 40 mm x 100 mm | 5x5x44" - 3.50 lb |
| ER51B | C-GL, POSPGL, Cuattro One | 1365 mm x 40 mm x 100 mm | 4x4x58" - 3.75 lb |
| ER51BT | 3 Panel Trans-Former | 1365 mm x 40 mm x 100 mm | 5x5x58" - 3.75 lb |
| ER34 | Any - Tape On | 960 mm x 40 mm x 3.2 mm | 4x4x42" - 2.50 lb |
| ER51 | Any - Tape On | 1315 mm x 40 mm x 3.2 mm | 4x4x56" - 2.75 lb |

FLAMMABILITY RATING

The product is non-flammable.
UL Rating: UL94V-2.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.
20148Milano Via don Gnocchi 24, Italy
+39 02 48706799

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:
Reina Imaging
6107 Lou Ave.
Crystal Lake, IL 60014 USA

DEVICE FAMILY CLASSIFICATION

| Country | Class |
|---------|--|
| USA | Class 1 per MDD 93/42/EEC, Annex V, NIST Compliant |
| Canada | Class 1 |
| Europe | Class 1, Annex VII, CE |
| TGA | Class 1 |

6107 West Lou Avenue | Crystal Lake, IL 60014 | (P) 800.752.4918 | (F) 815.356.8270 | reinaimaging.com