


## Product Datasheet

Model: PDRFDX3543  
Description: PROTECT-A-GRID DR FOR TOSHIBA FDX3543RP  
Revision: SEPTEMBER 2012



### Product Attributes:

- Holds DR Panel during exams.
- Encases grid to protect against damage.
- Provides additional protection for DR Panel.

### Dimensions

Product Dimension: 18.9 x 18.8 x 1 in. (48.1 x 47.8 x 2.54 cm)  
Weight: 2.2 lbs. (1 kg)

### Material Composition

Part	Material	Approximate %
Shell	HDP	98
Overlay	Lexan Film	2

### Weight Tolerance

*Note: This product was not designed for weight-bearing studies and should never be used without the DR Panel installed.*

Weight Tolerance: 400 lbs. (181.4 kg)

### Flammability Rating

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

### Biocompatibility

This model is a biocompatible device per ISO 10993.

### Statement of Conformance

Reina Imaging hereby states that the data, specifications, and descriptions contained herein in this document are accurate and representative of the models being described.

### Registrations and Compliance

All models listed above are made in the USA and registered to:

**Reina Imaging**

**6107 Lou Ave.**

**Crystal Lake, IL 60014 USA**

### Authorized EU Representative

PINCHET RX. S.r.l.

20148 MILANO, Via don Gnocchi 24, ITALY

### Device Family Classification

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1