


Product Datasheet

Model: PDRF560
Description: PROTECT-A-GRID DR FOR FLAATZ 560
Revision: MARCH 2012



Product Attributes:

- Holds DR Panel during exams.
- Encases grid to protect against damage.
- Provides additional protection for DR Panel.

Dimensions

Product Dimension: 18.9 x 18.8 x 1 in. (48.1 x 47.8 x 2.54 cm)
Weight: 2.2 lbs. (1 kg)

Material Composition

Part	Material	Approximate %
Shell	HDP	98
Overlay	Lexan Film	2

Weight Tolerance

Note: This product was not designed for weight-bearing studies and should never be used without the DR Panel installed.

Weight Tolerance: 400 lbs. (181.4 kg)

Flammability Rating

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

Biocompatibility

This model is a biocompatible device per ISO 10993.

Statement of Conformance

Reina Imaging hereby states that the data, specifications, and descriptions contained herein in this document are accurate and representative of the models being described.

Registrations and Compliance

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA

Authorized EU Representative

PINCHET RX. S.r.l.

20148 MILANO, Via don Gnocchi 24, ITALY

Device Family Classification

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1