


SIMPLY SOLID

HD Protect-A-Grid DRP Encasements HD Mobile DR Panel Holders



ReinaImaging


Reimagine Radiography

HD Protect-A-Grid DRP Encasements

PRODUCT ATTRIBUTES


- High tensile construction ideal for teaching hospitals, high volume facilities and bariatric patients
- Made from resilient, shock absorbent HDPE featuring radius edges for deflection and extra drop protection
- Curved corners and radius edges don't catch on bed sheets or clothes for easy placement
- Smooth, non-porous surface for effective cleaning and superior infection control
- Simple Twin-Lock™ dual locking design for easy panel application and removal
- Ergonomic handles to reduce fatigue and provide a sure grip
- Rugged one piece uniframe construction for durability
- Grid permanently protected inside sealed encasement

The **HD Protect-A-Grid DRP** grid encasement line provides additional DR panel protection for demanding applications. With its durable, heavy-duty composition, the HD line represents a powerful tool to protect expensive DR panels that are used under the most challenging conditions during free exposure exams. Key features include ergonomic handles for easy positioning in the most challenging environments, reliable Twin-Lock™ dual locking system and elegantly simple uniframe design for exceptional infection control and patient safety.



HD Tri-Handle Protect-A-Grid DRP grid encasement

55.2x45x2.54cm (21.75x17.75x1in.)
1.34kg (2.75lb. (w/out grid))



HD Dual Handle Protect-A-Grid DRP grid encasement

51.4x45x2.54cm (20.25x17.75x1in.)
1 kg (2.25lb. (w/out grid))

HD Mobile DR Panel Holders

PRODUCT ATTRIBUTES

- High tensile construction ideal for teaching hospitals, high volume facilities and bariatric patients
- Made from resilient, shock absorbent HDPE featuring radius edges for deflection and extra drop protection
- Curved corners and radius edges don't catch on bed sheets or clothes for easy placement
- Smooth, non-porous surface for effective cleaning and superior infection control
- Simple Twin-Lock™ dual locking design for easy panel application and removal
- Ergonomic handles to reduce fatigue and provide a sure grip
- One piece uniframe construction adds durability

The **HD Mobile DR Panel Holder** line provides additional protection and rigidity for demanding applications. With its durable featherweight composition, the HD line represents a powerful tool to protect expensive DR panels that are used under the most challenging conditions during free exposure exams. Key features include ergonomic handles configurable to user preferences for easy positioning in the most challenging environments, reliable Twin-Lock™ dual locking system and elegantly simple uniframe design for exceptional infection control and patient safety.



HD Tri-Handle Mobile DR Panel Holder

55.2x45x2.54cm (21.75x17.75x1in.)

1.34kg (2.75lb. (w/out grid))



HD Dual Handle Mobile DR Panel Holder

51.4x45x2.54cm (20.25x17.75x1in.)

1 kg (2.25lb. (w/out grid))

CONTACT REINA IMAGING

Every customer has unique requirements. Whether you are an OEM, dealer or end user, Reina Imaging has a dedicated team that can help. Together, we'll work to understand and create solutions for your particular needs. Give us a call or send us an email today to get started with your free consultation.

ADDRESS

6107 W. Lou Ave.
Crystal Lake, IL 60014

WEBSITE

reinaimaging.com

EMAIL

customer_service@reinaimaging.com

TOLL FREE

800-752-4918

PHONE

815-356-8181

FAX

815-356-8270



PAG
PROTECTA
2100 North Park Avenue, PA, CE

Ratio: 8:1
Focus: 4072"
LPI: 03
Date: 6/12/16
Model: P800
Serial: 1001200068734

Manufactured by Reina Imaging 9010 Lou St., Crystal Lake, IL 60014, USA
© 2016 Reina Imaging