

PRODUCT DATASHEET

Model: **P2D14-SS**

Description: **14X17 SURGICAL DROP-ON PROTECT-A-GRID 2 (ALUM) ENCASEMENT**

Revision: May 2018

ReinaImaging



PRODUCT ATTRIBUTES:

- Encases and protects grid.
- Fits over image receptor.
- Stainless steel 24" ruler with hook included for custom positioning under table.

MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Frame and Shims	Aluminum	98
Overlay	Polycarbonate Film	2
Ruler (separate part)	Stainless Steel	100

DIMENSIONS AND PERFORMANCE

Product Dimension: 18.5 x 15.0625 x 0.625 in. (47 x 38.26 x 1.6 cm)

Weight: 2.6 lbs. (1.18 kg) *without grid or ruler*

Weight Tolerance: 400 lbs. (181.4 kg)

Not designed for weight bearing exams. Not to be used without an image receptor.

FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and istered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

6107 West Lou Avenue | Crystal Lake, IL 60014 | (P) 800.752.4918 | (F) 815.356.8270 | reinaimaging.com