



Product Datasheet
Model: T34M-A, T34-A, T50M-A, T50-A
Revision: Original (July 2010)

Long-Leg / Full-Spine Grid Cabinet for use with Agfa CR or Film Cassettes

Product Attributes:

- Mobile or Wall-Mounted Installation
- Houses a removable grid cabinet for 2 or 3 Cassette Agfa Long-Leg / Full Spine Imaging Systems
- Grid encased in protective polycarbonate plastic to prevent incidental damage
- Cabinet adjusts vertically and tilts 90 degrees
- Non-electrical

Technical Specifications

Model	T34M-A (2 cassette)	T34-A (2 cassette)	T50M-A (3 cassette)	T50-A (3 cassette)
Weight (approx.)	88 lbs (39.9 kg)	88 lbs (39.9 kg)	86 lbs (39 kg)	86 lbs (39 kg)
Size	22.5 x 53.6 in.	22.5 x 53.6 in.	22.5 x 53.6 in.	22.5 x 53.6 in.
Vertical Travel (Height Adjustment)	Top - 74.25"	Top - 74.25"	Top - 75.75"	Top - 75.75"
Wall / Mobile	Mobile	Wall	Mobile	Wall
Boxed Dimension	31 x 31 x 60"	10 x 26 x 60"	31 x 31 x 60"	10 x 26 x 60"
Shipping Weight	120 lbs (54.4 kg)	75 lbs (34 kg)	120 lbs (54.4 kg)	75 lbs (34 kg)

Weight Tolerance

Note: All models above were not designed for performing standing feet x-ray exams and should never be used without cassettes installed.

Grid Cabinet Weight Tolerance = 600 lbs (272 kg) with cassettes installed

Flammability Rating

The product is non-flammable under the normal operating conditions. UL Rating: UL94V-2

Biocompatibility

All of the models listed in this data sheet are biocompatible devices.

Registrations and Compliance

Product registered to:

Reina Imaging
6107 Lou Ave.
Crystal Lake, IL 60014

- All models above are Class 1 Medical Devices in the USA
- CE
- 510(k) exempt
- Made in the USA

Contacts

Leo J Reina, President
815-356-8181
lreina@reinaimaging.com

Tony Reina, Quality Manager
815-356-8181
treina@reinaimaging.com