

## PRODUCT DATASHEET



Model: **PDRNG80C**

Description: **DIGITAL DETECTOR HOLDER FOR CANON CXDI-80C**

Revision: June 2017

### PRODUCT ATTRIBUTES:

- Holds DR Panel during exams.
- Protects DR Panel during routine handling.
- Provides ergonomic handle and double positive lock to secure DR Panel during transport.
- Made of shock absorbing materials.



### MATERIAL COMPOSITION

PART	MATERIAL	APPROXIMATE %
Shell	HDP	98
Overlay	Polycarbonate Film	2

### DIMENSIONS AND PERFORMANCE

Product Dimension: 15 x 16 x 1 in (38.1 x 40.64 x 2.54 cm)

Weight: 1.55 lbs. (.67 kg)

Weight Tolerance: 400 lbs. (181.4 kg)

*Not designed for weight bearing exams. Not to be used without an image receptor.*

### FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

### BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

### STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

### REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

**Reina Imaging**

**6107 Lou Ave.**

**Crystal Lake, IL 60014 USA**

### AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.

20148Milano Via don Gnocchi 24, Italy

+39 02 48706799

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1

6107 West Lou Avenue | Crystal Lake, IL 60014 | (P) 800.752.4918 | (F) 815.356.8270 | [reinaimaging.com](http://reinaimaging.com)