

Fully mechanical.  
**Motorized feel.**

# 80i Imaging Accessories

There are many challenges when radiographing long sections of anatomy. Patient hold times, motion, efficiencies, and outcomes are all impacted by the accessories holding or protecting the DR panel. Our full line of accessories designed to fit Fuji's 80i panel, match its sleek low-weight size, unlocking the continuous monolithic panel's true potential. Whether you are shooting an upright spine in the room or underneath a Jackson table intra-op, we have an accessory configuration that will achieve superlative results at the best price.



## 80i Mobile Positioning Partner

### PRODUCT ADVANTAGES

- Zero gravity lift assist
- Supine, landscape, upright orientations
- Easy to drape and store in surgery
- Telescoping horizontal arm
- 4 locking swivel casters
- Counter balanced weighted base
- Fully mechanical



Telescoping arm and locks



Supine Position



Landscape Profile

### Technical Data

Outside Dimensions	40.5 x 21 x 37.25 in. (102.87 x 53.34 x 94.62cm)
Total Vertical Stroke	43"/109.22cm
Max. Height Active Imaging Area - Portrait	75.5"/191.77cm
Min. Height Active Imaging Area - Portrait	1"/2.54cm
Max. Height Active Imaging Area - Landscape	23.25"/59.06cm
Min. Height Active Imaging Area - Landscape	6.5"/16.51cm
Weight (without 80I DR Panel and Grid mounted)	150lbs / 68Kg.



**Reina Imaging** CE

Safer • Protective • Preferred Radiology Accessories



[www.reinaimaging.com](http://www.reinaimaging.com)



[customer\\_service@reinaimaging.com](mailto:customer_service@reinaimaging.com)



## 80i Mobile Trans-Former

### PRODUCT ADVANTAGES

- Counter balanced vertical travel
- Rotates portrait and landscape
- Hospital grade casters
- Removable anti-scatter grid
- Fully mechanical, no electric needed (non-PBL)
- Wall Mounted version also available

### Technical Data

Dimensions (lowest portrait position)	53.5 x 18.5 x 23" (139.89cm x 47cm x 58.42cm)
Total Vertical Stroke	31.5"/80.01cm
Max. Height Active Imaging Area - Portrait	66"/167.64cm
Min. Height Active Imaging Area - Portrait	2.5"/6.35cm
Max. Height Active Imaging Area - Landscape	58.5"/148.59cm
Min. Height Active Imaging Area - Landscape	9.5"/24.13cm
Weight (without Long Length DR Panel and Grid mounted)	86lbs / 39Kg.



**ReinaImaging** CE


Safer • Protective • Preferred Radiology Accessories



[www.reinaimaging.com](http://www.reinaimaging.com)



[customer\\_service@reinaimaging.com](mailto:customer_service@reinaimaging.com)



## 80i Weight Distribution Cap

### PRODUCT ADVANTAGES

- Rated up to 500lb of distributed patient weight
- Scratch resistant and easy to clean polycarbonate vinyl overlay
- Slip resistant feet
- Cutouts for easy transport and positioning

## Technical Data

Outside Dimensions	35.1 x 19.7 x 1.5 in (89.2 x 50 x 3.81 cm)
Weight Capacity	500 lb. (227kg) distributed weight
Product Weight	14.35 lbs. (6.51 kg)

## 80i Removable Grid

### PRODUCT ADVANTAGES

- Scratch resistant polycarbonate overlay
- HDPE fully encased grid protection
- Includes highest quality 8:1/103lpi (40lpcm)/Parallel grid
- Other grid specifications available upon request



## Technical Data

Outside Dimensions	33.5 x 18.1 x .25 in (85.1 x 46 x .64 cm)
Product Weight	10.95 lbs. (4.97 kg)