

PRODUCT DATASHEET

Model: CC517, CC514, CC510, CC524

Description: CarbonLite™ Grid Encasements

Revision: June 2020



Product Attributes

- CarbonLite™ features .5mm carbon fiber
- Single ergonomic handle on long axis
- Rigid glass-filled polypropylene frame
- Curved corners and radiused edges
- Smooth, non-porous surfaces
- Simple Twin-Lock™ dual locking design
- Shock absorbent corner bumpers and inserts
- Fits all 17x17in, 14x17in, 10x12in and 24x30cm ANSI standard sizes

Material Composition

Material	%
Frame	
GF Polypropylene	85%
Aluminum	6%
Stainless Steel	5%
Flex PVC	2%
Acetal	1%
PC ABS	1%
Graphics Overlay	
Polycarbonate	98%
Colorant	2%
 Tubeside Protection	
Carbon Fiber	98%
Polycarbonate	2%

PRODUCT DATASHEET

Product Code	Product Name	Dimensions (cm)	Weight (kg)
17 x 17in			
CC517	17x17 CarbonLite Grid Encasement	56.42 x 49.22 x 2.69	1.76*
14 x 17in			
CC514	14x17 CarbonLite Grid Encasement	48.91 x 49.22 x 2.69	1.58*
10 x 12in			
CC510	10x12 CarbonLite Grid Encasement	38.47 x 36.54 x 2.69	1*
24 x 30cm			
CC524	24x30 CarbonLite Grid Encasement	38.47 x 36.54 x 2.69	1*

*Plus selected grid

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA