

CarbonMax

Product Advantages

- 2.5mm full carbon fiber tube side
- Single handle on long or short axis to help carry and position the image receptor
- Virtually indestructible, shock absorbent, glass-filled polypropylene frame
- Curved corners and radiused edges that won't catch on bed sheets or clothes for easy placement
- Smooth, non-porous surfaces for effective cleaning and superior infection control
- Simple Twin-Lock™ dual locking design for easy panel application and removal
- Ergonomic handle to reduce fatigue and provide a sure grip
- Shock absorbent corner bumpers and inserts for additional drop protection and perfect fit
- Fits most ANSI sized wireless DR panels and most tethered DR panels



The **CarbonMax™**, the toughest **CarbonCase** family member, features a full 2.5mm sheet of carbon fiber, providing the most substantial tube side protection. Available with and without a grid in a variety of sizes, the **CarbonMax** boasts a strength-to-weight ratio that makes it pound for pound the ***world champion in detector protection.***


Product Details




Tube Side



Rear Side with Twin-Lock™



Curved Corners and Cable Management



Shock Absorbent Corner Bumpers



ReinaImaging ce

Safer • Protective • Preferred Radiology Accessories



www.reinaimaging.com

 customer_service@reinaimaging.com

Product Specifications

Product Code	Product Name	Dimensions (cm)	Weight (kg)
17 x 17 Handle on Long Axis			
CC217	17x17 CarbonMax Grid Encasement	56.42 x 49.22 x 2.69	1.76*
CC2NG17	17x17 CarbonMax Mobile DR Panel Holder	56.42 x 49.22 x 2.69	1.76
14 x 17 Handle on Long Axis			
CC214	14x17 CarbonMax Grid Encasement	48.91 x 49.22 x 2.69	1.54*
CC2NG14	14x17 CarbonMax Mobile DR Panel Holder	48.91 x 49.22 x 2.69	1.54
10 x 12 Handle on Long Axis			
CC210	10x12 CarbonMax Grid Encasement	38.47 x 36.54 x 2.69	1*
CC2NG10	10x12 CarbonMax Mobile DR Panel Holder	38.47 x 36.54 x 2.69	1
24 x 30 Handle on Long axis			
CC224	24x30 CarbonMax Grid Encasement	38.47 x 36.54 x 2.69	1*
CC2NG24	24x30 CarbonMax Mobile DR Panel Holder	38.47 x 36.54 x 2.69	1

*Plus Selected Grid



ReinaImaging



www.reinaimaging.com



customer_service@reinaimaging.com

Safer • Protective • Preferred Radiology Accessories