

# Clear Advantage Stitching Stand

## Product Advantages

Height adjustable handles with easy-to-use squeeze locks

Four 360-degree swivel locking casters

Thin lightweight design

Radiolucent plexiglass support



Long leg and full spine imaging can be a challenge – especially in a digital environment. The **Clear Advantage Stitching Stand** is a critical tool for these specialized views. The simple and sleek design improves efficiencies and outcomes at an affordable price. The stand helps keep the patient immobile during multiple or single exposures. With four 360 degree lockable swivel casters, height adjustable handles and slim lightweight design, this immobilization device can be easily placed in position and removed for storage.

# Product Details



\*Trans-Former not included



**Reina Imaging** CE

Safer • Protective • Preferred Radiology Accessories



[www.reinaimaging.com](http://www.reinaimaging.com)

 [customer\\_service@reinaimaging.com](mailto:customer_service@reinaimaging.com)

# Product Details



\*Trans-Former and radiopaque ruler are not included

## Clear Advantage Stitching Stand

Outside Dimensions	30x29x79in. (76.2x73.6x200.7cm)
Handrail height (lowest)	21in. (53.3cm)
Handrail height (highest)	76in. (193cm)
Overall Weight	58lb. (26.31kg)



**Reina Imaging** CE

Safer • Protective • Preferred Radiology Accessories



[www.reinaimaging.com](http://www.reinaimaging.com)

 [customer\\_service@reinaimaging.com](mailto:customer_service@reinaimaging.com)