

DR IMAGING FOOTSTOOLS

CLEAR TWO STEP FOOTSTOOL

PRODUCT ADVANTAGES

- Ideal for lateral and AP weight bearing exams
- 3 separate slot positions for image receptor
- Safety handles for patient support
- Durable, non-slip surfaces
- 16 in. high platform
- Can accommodate a 17 x 17 in. cassette size imaging receptor or smaller
- Tilt-back casters for easy mobility
- Rated up to 500lb.



CLEAR SINGLE STEP FOOTSTOOL

PRODUCT ADVANTAGES

- Ideal for lateral and AP weight bearing exams
- 3 separate slot positions for image receptor
- Safety handles for patient support
- Durable, non-slip surfaces
- 9in. high platform
- Can accommodate a 17 x 17 in. cassette size imaging receptor or smaller
- Tilt-back casters for easy mobility
- Rated up to 500lb.

DR IMAGING FOOTSTOOLS



TWO STEP FOOTSTOOL

PRODUCT ADVANTAGES

- Ideal for lateral weight bearing exams
- Safety handles for patient support
- Durable, non-slip surfaces
- 16in. high platform
- size imaging receptor or smaller
- Tilt-back casters for easy mobility
- Rated up to 500lb

LOW BASE IMAGING FOOTSTOOL

PRODUCT ADVANTAGES

- Positions patient for weight bearing lateral foot and ankle exams
- Accommodates a 17x17in cassette size panel and smaller
- Slip resistant base for stability
- 9in. high platform
- Rated up to 500lb



Reina Imaging CE

Safer • Protective • Preferred Radiology Accessories



www.reinaimaging.com



customer_service@reinaimaging.com

DR IMAGING FOOTSTOOLS

U-ARM FOOTSTOOL PRODUCT ADVANTAGES

- Holds 500 lb.
- Standing feet exams
- Handrails for patient support
- Single step platform
- Optional height adjustable kit
- U-Arm Cavity: 10.25 x 28.5 in.



Reina Imaging CE

Safer • Protective • Preferred Radiology Accessories



www.reinaimaging.com



customer_service@reinaimaging.com