

EASY READER™ RADIOPAQUE RULER

Product Advantages

Effortless application

Lead-free, environmentally friendly,
safe for patient contact


Stainless steel plastic-lined bracket
slides easily across entire image area

Clearly visible, high-contrast, easy
to read radiopaque reference marks

Made from resilient material that is
both flexible and strong

Available in various lengths with and
without mounting bracket

The **Easy Reader radiopaque ruler** is ideal for spine and hips to ankle examinations. Its high-contrast NIST compliant reference marks are easy to see with your eyes as well as in a radiograph. It is simple to attach and remove with a stainless steel plastic-lined bracket that glides across the entire imaging area and won't scratch or mar equipment. Available in a variety of lengths, with or without a mounting bracket. It is 100% lead-free, RoHS compliant and made in the USA.



 **Rein Imaging** CE

Safer • Protective • Preferred Radiology Accessories




www.reinaimaging.com

 customer_service@reinaimaging.com

Product Details



Hook allows for easy application, positioning and removal



Floor mounted column with long format detector



Attaches to Stitching Platform