

Mobile DR Panel Full Enclosure

Product Advantages

IP 5X Rating

Shock absorbent foam lining

Quick release locks

Superior particulate and dust protection

Smooth non-porous surface

Curved corners

Available in 14x17in and 17x17in



The **14x17 Mobile DR Panel Full Enclosure** provides additional fluid and drop protection to sensitive flat panel detectors. With **quick release locks**, the case is both efficient and easy to use allowing the radiographer to focus on **positioning and patient care** instead of panel protection. Ideal for the ER, OR, disaster relief and any environment that is subject to potential fluid damage or accidental drops.

Product Details



Tube Side



Rear Side



Cutouts for LEDs and PB



Quick Release Locks

Technical Data for FE14

Product Dimensions

19.34 x 18.56 x 0.944 in. (49.12 x 47.14 x 2.39 cm)

Weight

3.20 lbs. (1.45 kg)



ReinaImaging ce

Safer • Protective • Preferred Radiology Accessories



www.reinaimaging.com



customer_service@reinaimaging.com