


HD Protect-A-Grid Encasements HD DR Panel Holders



The **HD Protect-A-Grid Encasement** and **HD DR Panel Holder** lines provide additional DR panel protection for demanding applications. With its durable, heavy-duty composition, the **HD line** represents a powerful tool to protect expensive DR panels that are used under the most challenging conditions during free exposure exams. Key features include ergonomic handles for easy positioning in the most challenging environments, reliable **Twin-Lock™ dual locking system** and elegantly simple uniframe design for exceptional infection control and patient safety.

HD Protect-A-Grid Encasements


Product Advantages

- High tensile construction ideal for teaching hospitals, high volume facilities and bariatric patients
- Made from resilient, shock absorbent HDPE featuring radius edges for deflection and extra drop protection
- Curved corners and radius edges don't catch on bed sheets or clothes for easy placement
- Smooth, non-porous surface for effective cleaning and superior infection control
- Simple Twin-Lock™ dual locking design for easy panel application and removal
- Ergonomic handles to reduce fatigue and provide a sure grip
- Rugged one piece uniframe construction for durability
- Grid permanently protected inside sealed encasement



**HD Dual Handle
Protect-A-Grid Encasement**
51.4x45x2.54cm (20.25x17.75x1in.)
1 kg (2.25lb. (w/out grid))

**HD Tri-Handle
Protect-A-Grid Encasement**
55.2x45x2.54cm (21.75x17.75x1in.)
1.34kg (2.75lb. (w/out grid))



HD DR Panel Holders

Product Advantages

- High tensile construction ideal for teaching hospitals, high volume facilities and bariatric patients
- Made from resilient, shock absorbent HDPE featuring radius edges for deflection and extra drop protection
- Curved corners and radius edges don't catch on bed sheets or clothes for easy placement
- Smooth, non-porous surface for effective cleaning and superior infection control
- Simple Twin-Lock™ dual locking design for easy panel application and removal
- Ergonomic handles to reduce fatigue and provide a sure grip
- One piece uniframe construction adds durability



HD Dual Handle Mobile DR Panel Holder

51.4x45x2.54cm (20.25x17.75x1in.)
1 kg (2.25lb. (w/out grid))

HD Tri-Handle Mobile DR Panel Holder

55.2x45x2.54cm (21.75x17.75x1in.)
1.34kg (2.75lb. (w/out grid))

