

PRODUCT DATASHEET

Model: **POSP14**

Description: **MOBILE DR PANEL POSITIONING PARTNER™**

Revision: March 2015



PRODUCT ATTRIBUTES:

- Telescoping arm
- Adjustable head
- Piston driven vertical travel
- Locking casters
- Counterbalanced

TECHNICAL SPECIFICATIONS

Weight (Approx.)	86 lbs. (39kg)
Size	30"x23.25"x40" (76.2cm x59cm x101.6cm)
Vertical Travel (Height Adjustment)	64" (162.56cm)
Horizontal Arm Vertical Travel	Max height from origin: 39" Min height from origin: 4"
Telescope Travel (Horizontal Arm Travel)	Max distance from origin: 29" Min distance from origin: 17"
Boxed Dimension	49" x 48" x 36" (124.5cm x 122cm x 91.5cm)
Shipping Weight	170 lbs. (86.2 kg)

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging
6107 Lou Ave.
Crystal Lake, IL 60014 USA

AUTHORIZED EU REPRESENTATIVE

PINCHET RX. S.r.l.
20148Milano Via don Gnocchi 24, Italy
+39 02 48706799

DEVICE FAMILY CLASSIFICATION

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1