

# PRODUCT DATASHEET

Model: POSPLC

Description: Positioning Partner Lite

Revision: January 2026

## Product Attributes

- Counterbalanced Vertical travel
- Small footprint
- Easy to drape
- No electric
- Locking casters
- Simple one-handed height adjustment

## Dimensions

Product Dimension: 28.3 x 18.8 x 58.3 in. (71.8 x 47.8 x 148.1 cm)

Total Vertical Stroke: 39 in. (99.06 cm)

Weight: 45 lbs. (20.4kg)

Boxed Dimensions: 24 x 20 x 57 in. (60.96 x 50.8 x 144.78 cm)

Boxed Weight: 75 lbs. (34kg)

## Flammability Rating

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

## Biocompatibility

This model is a biocompatible device per ISO 10993.

## Statement Of Conformance

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

## Registrations and Compliance

The model listed above is made in the USA and registered to:

Reina Imaging

6107 W Lou Ave.

Crystal Lake, IL 60014 USA

## Authorized Representative

Med Devices Lifesciences B.V.

Address: Keizersgracht 482, 1017 EG

Amsterdam, The Netherlands

Contact: +31-202254558

Email: info@meddevices.net

SRN: NL-AR-000000594



\*Protect-A-Grid Encasement not included