

PolyCarbon

Product Advantages

Single handle on long or short axis to help carry and position the image receptor

Virtually indestructible, shock absorbent, glass-filled polypropylene frame

Curved corners and radiused edges don't catch on bed sheets or clothes for easy placement

Smooth, non-porous surfaces for effective cleaning and superior infection control

Simple Twin-Lock™ dual locking design for easy panel application and removal

Ergonomic handle to reduce fatigue and provide a sure grip

Shock absorbent corner bumpers and inserts for additional drop protection and perfect fit

Fits most ANSI sized wireless DR panels and most tethered DR panels



The **PolyCarbon™ CarbonCase** provides a rugged, extremely durable polycarbonate top plate while offering the most affordable price in the market. Like all Reina Imaging detector holders, the **PolyCarbon** product line provides superior protection for both your DR panel and the grid resulting in the **lowest lifetime cost of ownership available today**.

Product Details




Tube Side



Rear Side with Twin-Lock™



Curved Corners and Cable Management



Shock Absorbent Corner Bumpers



Reinaimaging CE

Safer • Protective • Preferred Radiology Accessories



www.reinaimaging.com

 customer_service@reinaimaging.com

Product Specifications

Product Code	Product Name	Dimensions (cm)	Weight (kg)
17 x 17 Handle on Long Axis			
CC17	17x17 PolyCarbon Grid Encasement	56.42 x 49.22 x 2.69	1.85*
CCNG17	17x17 PolyCarbon Mobile DR Panel Holder	56.42 x 49.22 x 2.69	1.85
14 x 17 Handle on Long Axis			
CC14	14x17 PolyCarbon Grid Encasement	48.91 x 49.22 x 2.69	1.67*
CCNG14	14x17 PolyCarbon Mobile DR Panel Holder	48.91 x 49.22 x 2.69	1.67
10 x 12 Handle on Long Axis			
CC10	10x12 PolyCarbon Grid Encasement	38.47 x 36.54 x 2.69	1.08*
CCNG10	10x12 PolyCarbon Mobile DR Panel Holder	38.47 x 36.54 x 2.69	1.08
24 x 30 Handle on Long axis			
CC24	24x30 PolyCarbon Grid Encasement	38.47 x 36.54 x 2.69	1.08*
CCNG24	24x30 PolyCarbon Mobile DR Panel Holder	38.47 x 36.54 x 2.69	1.08

*Plus Selected Grid



Reina Imaging



www.reinaimaging.com



customer_service@reinaimaging.com

Safer • Protective • Preferred Radiology Accessories