

POP-UP CLINIC ACCESSORIES

Ideal diagnostic accessories for pop-up clinics in any environment. Our line of mobile and fixed positioning devices, image receptor holders and protective cases effectively address your general radiographic needs. With zero setup or installation and no requirement for power, our mechanical devices are service-ready out of the box and can withstand high-volume applications under the most adverse operating conditions.

MOBILE TRANS-FORMER® DRP HOLDER PRODUCT ADVANTAGES

- Ideal for upright Chest X-Rays
- No electric needed
- Zero gravity lift assist
- Can be configured to work with any 14x17in or 17x17in. image receptor
- Panel can be positioned portrait or landscape
- Vertical travel
- Pushbutton counterbalanced pulley system
- Removable grid
- AEC compatible
- Hospital grade casters



POP-UP CLINIC ACCESSORIES

MOBILE POSITIONING PARTNER®

PRODUCT ADVANTAGES

- Ideal for any upright or table top examinations
- Size adjustable head to hold 17x17in. image receptor or smaller
- Telescoping arm
- Adjustable head
- Piston driven vertical travel
- Small footprint
- Easy to drape
- No electric
- Locking casters
- Counterbalanced

View the video at [youtube.com/watch?v=xKeG-7OMM90](https://www.youtube.com/watch?v=xKeG-7OMM90)

The **Mobile DR Panel Positioning Partner** provides stability, versatility and reliability when positioning DR panels for portable radiography. Perpendicularity is crucial during cross-table exams and expensive panels are heavier than their imaging receptor predecessors - the Positioning Partner allows the tech to positioning effectively while achieving optimal images for a variety of radiographs. With its small footprint and counterbalanced base, the mobile DR panel holder is easy to drape for surgery and will not tip. It can be used in surgery, imaging centers, clinics, mobile pop-up clinics and exam rooms.



POP-UP CLINIC ACCESSORIES

PROTECT-A-GRID® DR PANEL ENCASEMENT

PRODUCT ADVANTAGES

- Made from resilient, shock absorbent HDPE featuring radius edges for deflection and extra drop protection
- Curved corners and radius edges don't catch on bed sheets or clothes for easy placement
- Smooth, non-porous surface for effective cleaning and superior infection control
- Simple Twin-Lock™ dual locking design for easy panel application and removal
- Ergonomic handles to reduce fatigue and provide a sure grip
- One piece uniframe construction adds durability



MOBILE DR PANEL FULL ENCLOSURE

PRODUCT ADVANTAGES

- Particulate and dust protection
- Drop protection
- Quick release locks
- Curved corners
- Smooth non-porous surface
- Lightweight
- Optional cutouts for LEDs and power button
- IP 5X

