

Positioning Partner Lite

Product Advantages

Counterbalanced Vertical travel

Small footprint

Easy to drape

No electric

Locking casters

Simple one-handed height adjustment



The **Positioning Partner Lite** features stability and reliability when positioning DR panels for cross table imaging. Its compact footprint makes it simple to drape and keep sterile when used in surgery and easy to store when not in use. The **Positioning Partner Lite** allows the technologist to position effectively and achieve optimal images for a variety of radiographs. It can be used in surgery, imaging centers, clinics and exam rooms. Built to the same rugged standards Reina Imaging is known for ensuring durability and long life.

*Protect-A-Grid Encasement not included

Product Details



*Protect-A-Grid Encasement not included

Whether you are positioning a cross table lateral c spine in the emergency department, a shoot through hip in the operating room or a lateral decubitus at the patient bedside, this is the way to ensure proper placement of the detector safely.

Easy one handed elevation, counterbalanced to ensure smooth operation.



Caster locks on all 4 wheels to secure the unit.



Product Details



*Protect-A-Grid Encasement not included