

# Stitching Platform With Handles

## Product Advantages

- Patient immobilizer during stitching exams
- Tilt-back casters for mobility
- Slip-resistant 9in. high platform
- Rated up to 500lb. (226.9kg) direct patient weight



LL and spine imaging can be tricky in a digital environment. The Stitching Platform with Handles is a critical tool for these specialized views. The simple design, keeps cost down while still performing as needed – keeping the patient immobile during multiple exposures and raising them up.



**Reina**Imaging 



Safer • Protective • Preferred Radiology Accessories



[www.reinaimaging.com](http://www.reinaimaging.com)



[customer\\_service@reinaimaging.com](mailto:customer_service@reinaimaging.com)



## Technical Data

Outside Dimensions	21 x 26.125 x 75.25 in. (53.34 x 66.36 x 191.14 cm)
Base Height	9 in. ( 22.86 cm)
Inside Width	24 in. (60.96 cm)
Overall Weight	61 lbs. (27.67 kg)
Weight Capacity	500 lb. (226.8 kg)
Biocompatibility	ISO 10993



**Reina Imaging** CE

Safer • Protective • Preferred Radiology Accessories



[www.reinaimaging.com](http://www.reinaimaging.com)



[customer\\_service@reinaimaging.com](mailto:customer_service@reinaimaging.com)