

# PRODUCT DATASHEET

Model: VFS

Description: DR Weight Bearing Foot Stool

Revision: November 2020

## Product Attributes

- For standing feet exams
- Handrails for patient support
- Single step platform
- Lightweight
- To be used with a U-Arm



## DIMENSIONS

Base Dimensions: 35x14x24 in. (88.9 x 35.56 x 60.96cm)

Handrail height: 35 in. (88.9 cm)

Weight: 33.6 lbs. (15.24 kg)

Weight Tolerance: 500 lbs. (226.8 kg)

U-Arm Cavity: 10.25 x 28.5 in. (26 x 72.4 cm)

Imaging Area: 17.25 x 17.25 in. (43.8 x 43.8 cm)

Boxed Dimensions and Weight

Box 1: 60 x 26 x 15" (152.4 x 66 x 38.1 cm), Boxed Weight: 75 lbs. (34 kg)

## FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

## BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

## STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

## REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA

| Country | Class                      |
|---------|----------------------------|
| USA     | FDA Class 1, 510(k) exempt |
| CANADA  | Class 1                    |
| EUROPE  | Class 1, Annex VII, CE     |
| TGA     | Class 1                    |